

Date:

Experiment:

Time fasted:

Time started:

Mouse	1	2	3	4	5	6	7	8	9	10	11	12	13	14	15
Genotype															
Age															
Sex															
Mass(g)															
Injection(μ L)															
Time(min) 0															

Notes:

Checklist:

Equipment for OGTT and ITT

Glucose meter

Glucose strips

Pen

Permanent marker

Calculator

Sharp scissor

Weight scale

Timer

20% glucose solution

(Vaseline)

For OGTT only

0.9% NaCl 0.5M EDTA in 96-well plate

Ice box

Heparinized capillaries

Feeding needle + 1mL syringes

200 μ L pipettor with tips

For ITT only

Insulin solution

Insulin syringes with needle for injection